

	Use Times (Shown In Hours)							
Flow Rate:	1	1.5	2	2.5	3	4	5	6
M4								
Pulse Dose	5.7	3.8	2.9	2.3	1.9	1.4	1.1	.9
Continuous Flow	1.9	1.3	.9	.7	.6	.5	.4	.3
M6								
Pulse Dose	8.3	5.5	4.1	3.3	2.8	2.1	1.7	1.4
Continuous Flow	2.7	1.8	1.4	1.1	.9	.7	.6	.4
ML6								
Pulse Dose	8.6	5.7	4.3	3.4	2.9	2.1	1.7	1.4
Continuous Flow	2.8	1.9	1.4	1.1	.9	.7	.6	.4
C								
Pulse Dose	12.1	8.1	6.1	4.9	4.0	3.0	2.4	2.0
Continuous Flow	4.0	2.7	2.0	1.6	1.3	1.0	.8	.7
D								
Pulse Dose	21.0	14.0	10.5	8.4	7.0	5.2	4.2	3.5
Continuous Flow	6.9	4.6	3.5	2.8	2.3	1.7	1.4	1.2
E								
Pulse Dose	34.4	23.0	17.2	13.8	11.5	8.6	6.9	5.8
Continuous Flow	11.4	7.6	5.7	4.6	3.8	2.8	2.3	1.9